



MEDICAL REVITALIZATION, INC

Hill Rom Transtar P8010



Features

- New Mattress
- New Paint
- Pumps and Stickers

Refurbishment Process

- Disassemble Bed
- Steam clean frame and components
- Paint frame
- Install new brake mechanisms on both sides
- Sand and repaint shrouds with our special 2 part polyurethane coating

Specifications

- Height Overall: High Position: 34", Low Position 21.50"
- Overall Length: 92"
- Overall Width: 26.50" (Siderails Up 32")
- Overall Siderail Length: 47"
- Caster Size: 6"
- Maximum Back Incline: 90°
- Maximum Weight Limit: 500lbs

- Rebuild pumps with new seals and adjust pumps speed and travel
- New stickers per customers choice of color
- Clean,adjust and lube brake/steer components
- New bumpers and plastic
- 1 year warranty