



MEDICAL REVITALIZATION, INC

Stryker 1001



Features

-

Specifications

Refurbishment Process

- Disassemble Bed
- Steam clean frame and components
- Paint frame
- Install new brake mechanisms on both sides
- Sand and repaint shrouds with our special 2 part polyurethane coating
- Rebuild pumps with new seals and adjust pumps speed and travel

- Overall Length: 83"
- Overall Width: Siderails Up 31.50", Siderails Down 30"
- Overall Height: High Position 35.50", Low position 22"
- Patient Surface: 26" x 76"
- Siderails: 14" x 54"
- Caster Size: 8"
- Knee Gatch: 0° to 30°
- Backrest: 0° to 90°

- New stickers per customers choice of color
- Clean,adjust and lube brake/steer components
- New bumpers and plastic
- 1 year warranty

- Trendelenburg/Reverse Trendelenburg: $\pm 18^\circ$

- Maximum Patient Weight: 500lbs