



MEDICAL REVITALIZATION, INC

Hill-Rom Procedural Recliner



Features:

- A manually adjustable back section which is controlled by a spring assisted release handle, located on both the left and right sides of the table.
 - The back section is adjustable from a full horizontal table position to a chair position.
- An extendable footrest shelf which includes a positionable upholstered pad and a removable treatment pan.
- Right and left stirrups which can be extended to an infinite number of positions as well as laterally in four positions.
- Two storage drawers on steel ball bearing glides which are located on the foot end of the table and a bulk storage area accessed through two hinged doors on both sides of the table.
 - An optional pelvic lift bar which allows the seat section to be positioned in the Lithotomy position.
- A self-locking steel step with a non-skid surface. The steps extends from the foot end of the table to assist patient entry and exit.
 - A duplex, hospital grade 115 VAC receptacle, located on the foot end, left side of the table.
 - An optional drawer warmer which warms the foot end, top drawer.
- A paper roll holder which is hidden in the head end of the upholstered top and tear strips, located at the foot end and head end of the upholstered top.