



MEDICAL REVITALIZATION, INC

Midmark 622



Features

- Heated seat and lumbar to improve patient comfort and reduce anxiety
- Upholstery and drawer heated controls are easily accessible for you and your staff
- Unique elastomeric material stretched over a seat cut-out is designed to enhance patient comfort for a "chair-like feel" in an exam table
- Allows for color change in seconds or additional access for a thorough cleaning. No tools needed.
- The Barrier-Free™ examination tables have the option of two easy-to-use foot controls, which

Specifications

- Top Dimensions: 28.38"W x 59.50"L
- Length with footrest extended: 77.50"
- Height: Minimum 18", Maximum: 37"
- Patient Load Rating: 400lbs.
- Paper Roll Holder: 21" x 3.50" Paper Roll
- Electrical Receptacle: 115VAC
- Rear Storage Pod: 17.80" x 7.0" x 9.80" High (includes removable divider)

power your patient effortlessly to the desired position. The 622 foot control raises and lowers the table to your requirements.

Refurbishment Process

- Disassemble Chair
- Remove Motors
- Steam clean internal Parts
- Replace gas assisted springs
- Computer test motors and circuit boards
- Lube gears
- Repaint all
- Polish plastics
- Recover fabric
- 1 Year Warranty