



MEDICAL REVITALIZATION, INC

Hill Rom Total Care ICU Bed



Features

- Continuous Lateral Rotation Therapy
- Percussion and Vibration Therapy, including postural option
- Weight based pressure redistribution in any bed position
- Advanced Microclimate Management technology
- Point-Of-Care® siderail controls with Graphical Caregiver Interface feature (GCI®)
- Digital Head of bed angle indicator and alarm
- 3 mode bed exit sensor with alarm silence and

Specifications

- Overall Length: FlexAfoot Extended 93", FlexAfoot Retracted 81.5"
- Overall Width: Siderails Stowed 38", Rails Up 40"
- Max Head Elevation: 75°
- Maximum Weight Capacity: 500lbs

alarm suspend modes

- FullChair® and chair egress with stand assist
- One-button Boost® patient repositioning support
- Turn assist
- Optional IntelliDrive® powered transport
- History data on patient weight, head of bed elevation, chair position, and pulmonary therapies

RefurbishmentProcess

- Disassemble Bed
- Steam clean frame and components
- Paint frame
- Install new brake mechanisms on both sides
- Sand and repaint shrouds with our special 2 part polyurethane coating
- Rebuild pumps with new seals and adjust pumps speed and travel
- New stickers per customers choice of color
- Clean,adjust and lube brake/steer components
- New bumpers and plastic
- 1 year warranty