



MEDICAL REVITALIZATION, INC

Hill Rom Versa Care



Features

- Boost one-button feature helps patients move toward the head of the bed
- IntelliDrive powered transport lets staff handle the bed effortlessly
- Point of Care control panel uses easy to read digital display and simple icons
- Fully featured low bed can be adjusted to the optimal height for patient care with full functionality in the lowest position
- Hands free foot pedal system allows caregivers to raise or lower the head of the bed and height of

Specifications

- Bed Height: High Position 34", Low Position 17"
- Overall Width: 40" (Siderails Down: 37")
- Overall Length: Fully Extended 94.50", Fully Retracted 82.50"
- Max Head Elevation: 0° to 65°
- Max Foot Section: 0° to 27°
- Max Patient Weight Limit: 500lbs

the bed without using their hands

- Low chair position allows patients to be placed in an upright position while keeping the bed low

Refurbishment Process

- Disassemble Bed
- Steam clean frame and components
- Paint frame
- Install new brake mechanisms on both sides
- Sand and repaint shrouds with our special 2 part polyurethane coating
- Rebuild pumps with new seals and adjust pumps speed and travel
- New stickers per customers choice of color
- Clean, adjust and lube brake/steer components
- New bumpers and plastic
- 1 year warranty